
Intended Use

IL Test™ AST/GOT is intended for the quantitative *in vitro* diagnostic determination of aspartate aminotransferase (AST) or glutamic oxalacetat transaminase, GOT) in human serum and plasma using the ILab Chemistry Systems. AST measurements are used in the diagnosis and treatment of certain types of liver and heart diseases. Refer to the publication "Effects of Diseases on Clinical Laboratory Tests;" for a summary of causes of increased and decreased AST.

Principle/Test Summary

Rate analysis. The reagent formulation is according to the recommendations of the IFCC.²

The rate of decrease in absorbance due to the oxidation of NADH to NAD⁺ is directly proportional to sample AST activity. Absorbance measurements are taken at 340 nm. The lactate dehydrogenase (LDH) included in the reagent eliminates interference from serum pyruvate.

Reactive Ingredients of Final Assay Mixture

| | |
|--|----------------|
| - 2-oxoglutarate | 12,0 mmol/L |
| - L-aspartate | 240 mmol/L |
| - Reduced nicotinamide adenine dinucleotide (NADH) | > 0,180 mmol/L |
| - Lactate dehydrogenase (LDH) microbial | > 600 U/L |
| - Malate dehydrogenase (MDH) microbial | > 420 U/L |
| - Buffer Tris | pH 7.7 |

Precautions

Contains sodium azide (< 0.1%) that may form explosive azides in metal plumbing. Dilute with water prior to disposal. During disposal into sewerage flush with copious amounts of running water. Metal plumbing should periodically be decontaminated with 10% sodium hydroxide. Handle all specimens as potentially infectious for HIV or Hepatitis B. Avoid contact with skin and eyes (S24/25). Do not empty into drains (S29). Wear suitable protective clothing (S36). This product is For *in vitro* Diagnostic Use.

Reagent preparation

AST/GOT R1: Reconstitute the AST/GOT R1 container with 20 mL of AST/GOT Buffer. Do not shake or allow to foam. Allow the working reagent to stand for 30 minutes before use. Place into reagent tray (ILab 900/1800; place into reagent tray #1). Note: IFCC reference method recommends the use of pyridoxalphosphate as activator. Pyridoxalphosphate auxiliary reagent (PLP) is available separately (Cat. No. 0018162750). If this method is used, add 1.0 mL of prepared PLP to the working AST/GOT R1 reagent. AST/GOT R2: Reagent is ready-to-use. Place in reagent tray (ILab 900/1800; place in reagent tray #2).

Storage of Unopened Reagents

Refrigerated (2-8°C) until expire. See expiration date on the carton.

Working Reagent Storage and Stability

10 days on board with pyridoxalphosphate.
15 days on board without pyridoxalphosphate.

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Indicators of Reagent Deterioration

Physical appearance: All reagents should be clear and colorless. *Control Assay:* The proper range should be obtained during a controlled assay. *Absorbance Range:* Changes in typical reagent blank absorbances may indicate reagent deterioration.

Specimens

Serum, lithium heparinized plasma. Samples should be analyzed as soon as possible after collection or stored refrigerated (2-8°C) for up to 24 hours.² Samples should be brought to room temperature and mixed before analysis.

Calibration

Enzyme activity is measured by the rate of production of a chromophore with well defined optical properties. The molar absorptivity constant of NADH is 6.3 x 10³. Maintenance of optical properties, fluidics and thermal control is critical to the accuracy of this assay. The ILab continuously monitors fluidics, thermal control and some optical parameters. Automatic daily maintenance monitors additional optical properties. Maintenance procedures are detailed in the ILab Operator's manual. Any adjustments to the ILab should be immediately followed by performance of daily maintenance and quality control procedures. Any failures in the continuous monitoring, daily maintenance or quality control procedures indicate improper system performance. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

Data Analysis

The ILab calculates AST activity as follows:

$$ACT_{sample} (U/L) = EF \times \Delta A \quad EF = \frac{VR}{V_s} \times \frac{1}{E} \times \frac{1}{PL} \times 10^6$$

Where
ACT = activity
EF = Enzyme factor
V_s = Total sample volume (mL)
PL = Pathlength of cuvette (cm)
* The measurement in the ILab cuvette is (0.5 cm for ILab 600 and 0.6 cm for ILab 900/1800) automatically converted to 1 cm pathlength by the ILab.

For AST/GOT in the ILab 900/1800:

$$EF = \frac{0.305}{0.025} \times \frac{1}{6.3 \times 10^3} \times \frac{1}{1} \times 10^6 = 1936$$

For AST/GOT in the ILab 600:

$$EF = \frac{0.194}{0.014} \times \frac{1}{6.3 \times 10^3} \times \frac{1}{1} \times 10^6 = 2200$$

Activity is expressed in International Units/Liter (U/L). To convert U/L to nKat/L, multiply by 16.67.³

Quality Control

Use normal and abnormal controls in accordance with good laboratory practices. Analyze at least once per day. Use at least 2 levels of IL SeraChem® controls, such as SeraChem® Level 1 (Cat. No. 0018162412) and SeraChem® Level 2 (Cat. No. 0018162512). The target values and expected range for this control materials can be found in the insert sheets. For identification and resolution of out-of-control situations, references such as Westgard *et al.*⁴ are recommended.

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Interferences

ILab 600: No interference from lipemia up to sample absorbance of 5.0/cm at 660 nm (970 mg/dL triglycerides). No interference up to 70 mg/dL (0,042 mmol/L) hemoglobin. No interference up to 30 mg/dL (510 µmol/L) bilirubin. ILab 900/1800: No interference from lipemia up to sample absorbance of 1,8/cm at 660 nm. No interference up to 125 mg/dL (0,075 mmol/L) hemoglobin. No interference up to 20 mg/dL (340 µmol/L) bilirubin. For a comprehensive review of interfering substances, refer to the publication by Young.⁵

Reference Range

Males: 15 to 48 U/L (250 - 800 nKat/L).⁶ Females: 9 to 36 U/L (150 - 600 nKat/L).⁶ These reference ranges assume the use of PLP. IL recommends that each laboratory establish its own reference range since ranges may vary with age, diet, gender, and geographical area.

Precision

The coefficient of variation (3 samples/run, 20 runs) for a sample with 20 U/L AST should not be greater than 5% within run and 7% total. Typical precision results obtained according to NCCLS protocol:⁷

| ILab 600 | Samples/Runs | Mean (U/L) | CV (%) | Mean (U/L) | CV (%) |
|---------------|--------------|------------|--------|------------|--------|
| Within Run | 3/20 | 25 | 2.3 | 52 | 1.0 |
| Total | 3/20 | 25 | 3.8 | 52 | 1.2 |
| ILab 900/1800 | | | | | |
| Within Run | 3/20 | 13 | 3.8 | 50 | 1.3 |
| Total | 3/20 | 13 | 4.3 | 50 | 1.5 |

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- Friedman R., Young D.S., Effects of Disease on Clinical Laboratory Tests, AACCPress, Washington, D.C., 1989.
- Bergmeyer H.U., *et al.*, "Approved Recommendations on IFCC Methods for the Measurements of Catalytic Concentrations of Enzymes, Part 2, IFCC Method for Aspartate Aminotransferase", *J. Clin. Chem. Clin. Biochem.*, 24:481, 1986.
- Lippert H., and Lehmann H.P., SI Units in Medicine, Urban and Schwarzenberg, Baltimore, 1978, p.86.
- Westgard J. and Barry P., Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACCPress, Washington, D.C., 1986.
- Young D.S., Effects of Drugs on Clinical Laboratory Tests, AACCPress, Washington, D.C., 1990.
- Tietz N., Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia, PA., 1990, pp 78.
- Document EP3-T: "Guidelines for Manufacturers for Establishing Performance Claims for Clinical Chemical Methods, Replication Experiment Evaluation", National Committee for Clinical Laboratory Standards, Villanova, PA, 1982.

iL Test™, SeraChem® are trademarks of Instrumentation Laboratory Company.

IVD

In Vitro Diagnostic Medical Device
Dispositivo medico-diagnostico in vitro
In Vitro Diagnostikum
Producto sanitario para diagnóstico in vitro
Dispositif medical de diagnostic- que in vitro
Προϊόν για διαγνωστική χρήση *in vitro*

LOT

Batch code
Codice del lotto
Chargenbezeichnung
Codigo de lote
Code du lot
Αρ. Παρτίδας

Use by
Utilizzare entro
Verwendbar bis
Fecha de caducidad
Utiliser jusque
Χρήση έως

Temperature limitation
Limite di temperatura
Zulässiger Temperaturebereich
Limite de temperatura
Limites de température
Περιορισμοί θερμοκρασίας

Consult instructions for use
Consultare le istruzioni per l'uso
Gebrauchsanweisung beachten
Consulte las instrucciones de uso
Consulter les instructions d'utilisation
Συμβουλευτείτε τις οδηγίες χρήσης

Manufacturer
Fabbricante
Hersteller
Fabricante
Fabricant
Κατασκευαστής

| iLab 600 Test Parameters | | |
|-----------------------------------|-------------------------|---------------------|
| Photometric Test Parameters | | Serum |
| Test No. | | 31 |
| Test Name, Test Code | | AST/GOT, AST |
| Sample Type | | Serum |
| Reporting Unit, Decimal Points | | U/L (nKat/L), 0 (0) |
| Reaction Cycle | | Standard |
| Twin Analysis | | OFF |
| Methodology Type, Measuring Point | | Rate, 20/27 |
| Photometric Methodology | | 1 Wavelength |
| Primary/Secondary Wavelength | | 340 |
| Sampling Conditions | | |
| Sampling 1 | Sample Vol. | 14 |
| | Sample/Diluent Vol. | 0/0 |
| Sampling 2 | Sample Vol. | 7 |
| | Sample/Diluent Vol. | 0/0 |
| Sampling 3 | Sample Vol. | 14 |
| | Sample/Diluent Vol. | 10/90 |
| Sampling 4 | | *** |
| Diluent Code | | Water |
| Diluent Warning Limit | | N/A |
| First Run | | Sampling 1 |
| Below/Above Normal Range | | *** |
| Panic L | | *** |
| Panic H | | Sampling 2 |
| Noise | | *** |
| Prozone | | N/A |
| HIGH!, ABS! | | Sampling 3 |
| Reagent Volumes | | |
| R1 | Code | 01311 |
| | Rgt/Dil. Vol., Stirring | 70/75, ON |
| | Low Vol. Warning Limit | *** |
| | Stability (days) | 10 |
| R2 | Code | 01312 |
| | Rgt/Dil. Vol., Stirring | 35/0, ON |
| | Low Vol. Warning Limit | *** |
| | Stability (days) | 10 |

| Ranges and Evaluation Criteria | Serum |
|---------------------------------|--------------------------------|
| Normal Range-Male | 15 - 48 (250 - 800) |
| Normal Range-Female | 9 - 36 (150 - 600) |
| Normal Range-Other | ** |
| Valid Range | 3 - 400 (50 - 6668) |
| Hemolysis/Icterus/Lipemia Limit | *** |
| Reaction Slope | Negative |
| Absorbance Limit | Below, 1000 |
| Prozone Limit | N/A |
| Non Linear Limit | 50 |
| Slope/Intercept Correction | 1/0 |
| Qualitative Report | OFF |
| Calibration Conditions | |
| Calibration | Factor: -2.200 (-36.7), 3 Reps |
| Stability (days) | N/A |
| Calibrator, Concentration | N/A |
| R-Blank Limit (mAbs) | 0 |
| Cal. Reps Range (%) | N/A |
| Min. Cal. Response (mAbs) | N/A |
| Cal. Factor Change (%) | N/A |
| M-Point Curve Fit (%) | N/A |
| Reagent Blank | ON |
| Auto R-Blank by Bottle | ON |

* Lot dependent - see package insert
** Operator definable
***Optional
N/A Not applicable to this test

| iLab 900/1800 Test Parameters | | | | |
|-------------------------------|-----------------|---------------|-------------|-------|
| Test Name | [AST] | | | |
| Bar code | [0131] | | | |
| Measuring method | [Rate] | | | |
| Calibration method | [Inst.Const.] | | | |
| Unit | [U/L] nKat/L | | | |
| Measuring points | [20] | [30] | [] | [] |
| Photometric method | [Single] | | | |
| Wavelength 1,2,3 | [340] | [] | [] | |
| Sample Vol. N, D, I | [25] | [10] | [30] | |
| R1V-Stir Mode [Dil.] | [110] | [Yes] | [120] | |
| R2V-Stir Mode [Nor.] | [50] | [Yes] | [] | |
| No. decimal points | [0] | | | |
| Sample type | [Serum] | | | |
| Normal range | MALE L-H | [15] 250 | [48] 800 | |
| | FEMALE L-H | [9] 150 | [36] 600 | |
| Valid range LL-HH | [3] 50 | [400] 6668 | | |
| Correction const. 1, 2 | [1.0] | [0.0] | | |
| W1-Abs Inc/Dec, Limit | [Dec] | [1000] | | |
| Hemolysis limit | [+] | | | |
| Icteric limit | [+] | | | |
| Lipemic limit | [+] | | | |
| Calib. Replicates | [4] | | | |
| Replicates, Sens. Limit *** | [2500] | [0.0] | | |
| Slope-Offset limits *** | [999] | [2500] | | |
| Partner Test, Const. | [++] | [++] | | |
| Prozone Judge, Limit | [++] | [++] | | |

K1 (Slope) constant [-1.936] -32.3

*** Operator definable data
+ See the interference section of the insert sheet
++ Not applicable to this method

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